

GE diffusion gradient directions (CV11: custom tensor file, tensorXX.dat)

Jaemin Shin (Jaemin.Shin@ge.com)

Updated: May 13, 2021



Imagination at work.

Multi-shell DTI or DSI: How to prescribe?

With a Max b-value, multi-shell DTI (or any smaller b-values) can be specified using the entries (b-vector, mag<1) of custom tensor filename (CV11).

$$\text{b-value}_{\text{effective}} = \text{b-value}_{\text{max}} * (\text{b-vecX}^2 + \text{b-vecY}^2 + \text{b-vecZ}^2)$$

2000 s/mm²


$$\begin{aligned} \text{b-value}_{\text{effective}} &= 2000 \times 0.387298^2 \\ &= 300 \text{ s/mm}^2 \end{aligned}$$

$$\begin{aligned} \text{b-value}_{\text{effective}} &= 0 \text{ s/mm}^2 \end{aligned}$$

	b-vecX	b-vecY	b-vecZ
107	0.387298	0.000000	0.000000
	0.000000	0.387298	0.000000
	0.000000	0.000000	0.387298
	-0.272813	-0.034121	-0.272813
	0.272813	0.272813	0.034121
	0.034121	0.272813	0.272813
	0.000000	0.000000	0.000000
	-0.632456	0.000000	0.000000
	0.000000	0.632456	0.000000
	0.000000	0.000000	0.632456
	0.026816	-0.072479	-0.627712
	-0.110616	-0.126807	-0.609687
	-0.146919	-0.102837	-0.606525

/usr/g/bin/tensor107.dat

b-value
definition



$$b = \gamma^2 G^2 \delta^2 (\Delta - \delta/3)$$

Details Diffusion Advanced

b-values: 1

#	b-value	NEX
1	2000.0	1.0

NEX for T2: 1.0

Diffusion Direction: TENSOR

of Diffusion Directions: 107

of T2 Images: 1

Synthetic b-values: 0

Details Diffusion Advanced

User Control Variables

		Min	Max
CV5	Recon Type (0=Zero Filling, 1=Homodyne):	1.00	0.0 1.0
CV9	Shim Volume Mode (0=Default, 1=Breast):	0.00	0.0 2.0
CV11	Tensor filename = tensor[n].dat (research use >99):	107.00	0.0 65536.0



GE convention for diffusion gradient direction

- GE convention for diffusion gradient direction has always been in ***“MR physics” logical coordinate***
 - X: Freq
 - Y: Phase
 - Z: Sliceencoding direction
- This is **neither** "with reference to the scanner bore" (like Siemens or Philips) **nor** "with reference to the imaging plane" (as expected by FSL tools)

tensor.dat

```
60
1.000000 0.000000 0.000000
0.361000 0.933000 0.000000
-0.255000 0.565000 0.785000
0.861000 -0.464000 0.210000
-0.307000 -0.766000 0.564000
-0.736000 0.013000 0.677000
0.532000 0.343000 0.774000
0.177000 0.965000 0.195000
0.771000 0.163000 0.615000
```

X: Freq Y: Phase Z: Slice

(0019,10bb) (0019,10bc) (0019,10bd)

DICOM private tags

dc2niix



FSL bval/bvec format
(with reference to the imaging plane)



space-separated vs tab-separated values in tensorXX.dat

- Software version prior to **29.1**,
- PSD: played-out diffusion gradient is applied **correctly** regardless of “spaces” vs “tab” separator.
- **However**, information in DICOM tags is **lost (wrong)** with “tab” separator. (X=Y=Z in DICOM)
- Version 29.1 or later, you can use “tab”

space-separated

```
102
0.000000 0.000000 0.000000
0.654875 0.355659 0.666817
0.271924 0.933965 0.231877
-0.781703 0.152743 -0.179657
-0.118810 0.110437 0.986756
-0.552748 0.108005 -0.127037
-0.326879 0.866547 0.377155
0.394817 0.108337 -0.912350
0.155639 0.477851 0.643508
-0.906253 0.349334 -0.238058
0.110053 0.337891 0.455029
-0.566330 0.470917 -0.676393
0.558540 -0.665203 -0.495518
0.002542 0.026310 -0.407392
-0.648904 -0.611989 -0.452098
0.265763 -0.689599 -0.347116
-0.435216 0.497941 0.750095
0.060085 -0.972644 -0.224396
-0.187923 0.487620 0.245448
```

tensor4321.dat

X: (0019,10BB): 0.654875
Y: (0019,10BC): 0.355659
Z: (0019,10BD): 0.666817

tab-separated

```
102
0.000000 0.000000 0.000000
0.654875 0.355659 0.666817
0.271924 0.933965 0.231877
-0.781703 0.152743 -0.179657
-0.118810 0.110437 0.986756
-0.552748 0.108005 -0.127037
-0.326879 0.866547 0.377155
0.394817 0.108337 -0.912350
0.155639 0.477851 0.643508
-0.906253 0.349334 -0.238058
0.110053 0.337891 0.455029
-0.566330 0.470917 -0.676393
0.558540 -0.665203 -0.495518
0.002542 0.026310 -0.407392
-0.648904 -0.611989 -0.452098
0.265763 -0.689599 -0.347116
-0.435216 0.497941 0.750095
0.060085 -0.972644 -0.224396
-0.187923 0.487620 0.245448
```

tensor4321.dat

X: (0019,10BB): 0.654875
Y: (0019,10BC): 0.654875
Z: (0019,10BD): 0.654875

Wrong



Information in DICOM & Software version

Diffusion b-value:

Public tag (0018, 9087), Private tag (0043, 1039)

Diffusion b-vector:

X: (0019, 10bb), Y: (0019, 10bc), Z: (0019, 10bd)

For software version 26 or older (ex. DV26),

- If you modify PSD to use a custom tensor file (for example, tensor1234.dat), PSD played-out diffusion gradient is correctly from tensor1234.dat.
- However, information in DICOM (0019, 10bb/c/d) is from tensor.dat, **not** from your custom tensor1234.dat specified in PSD

For software version 27 or later,

- Both PSD played-out diffusion gradient and information in DICOM is from the custom filename (CV11)




Multi-shell DTI in clinical mode (version 29.1)

- From software version **29.1**, custom tensor (CV11) in clinical mode
- Pre-load most widely used tensor files in product
 - ADNI, HCP, ABCD, UK biobank, ABCD
- Correct b-value in DICOM tags (0018, 9087) (0043, 1039)
 - Round to nearest multiple of 5 (ex. 1999.8 → 2000)
- dcm2niix** (v1.0.20201102 or later) reports bval/bvec correctly and consistently for both new/existing datasets

Correct b-value
annotation (ABCD)

b-value
[s/mm²]



Details Diffusion Advanced				
User Control Variables				
			Min	Max
CV5	Recon Type (0=Zero Filling, 1=Homodyne):	1.00	0.0	1.0
CV9	Shim Volume Mode (0=Default, 1=Breast):	0.00	0.0	2.0
CV11	Tensor filename = tensor[n].dat:	90.00	0.0	65536.0



Adolescent Brain Cognitive Development®
Teen Brains. Today's Science. Brighter Future.

Appendix:

Pre-loaded multi-shell tensor files



Imagination at work.

tensor8313: HCP (Human Connectome Project)

tensor8313.dat

- Three-shells DTI
- # of directions = 270, Max b-value = 3000
 - 1000 x 90
 - 2000 x 90
 - 3000 x 90
- # of directions = 296, Max b-value = 3000
 - WITH interleaved b0 values
 - 0 x 26
 - 1000 x 90
 - 2000 x 90
 - 3000 x 90




```
# Multi-shell tensor8313.dat file
#
# Date generated: 2018-02-07 14:33:01
#
# Please use the following configuration on the MR Operator Console:
# - Advanced tab: CV11 - tensor file name: n=8313
# - Diffusion tab: b-value: 3000
#
# - WITH interleaved b0 values (270 directions + 26 b0 volumes):
#   - Diffusion tab: number of directions (TENSOR): 296
#
# - WITHOUT interleaved b0 values:
#   - Diffusion tab: number of directions (TENSOR): 270
#
#
# This file should be placed in /usr/g/bin
#
# Multi-shell details:
# - DTI shell 1: b value = 1000; directions = 90
# - DTI shell 2: b value = 2000; directions = 90
# - DTI shell 3: b value = 3000; directions = 90
#
# Interleaved b0 details:
# - The 296 directions comprise 270 tensor directions and 26 b0 volumes, each one interleaved every 10 directions.
#
# Identical directions are unique.
```



tensor4321: ABCD (Adolescent Brain Cognitive Development)

tensor4321.dat

- Four-shells DTI
- # of directions = 96, Max b-value = 3000
 - 500 x 6
 - 1000 x 15
 - 2000 x 15
 - 3000 x 60
- # of directions = 102, Max b-value = 3000
 - WITH interleaved b0 values
 - 0 x 6
 - 500 x 6
 - 1000 x 15
 - 2000 x 15
 - 3000 x 60



```
# Multi-shell tensor4321.dat file for ABCD (Adolescent Brain Cognitive Development) study
# Diffusion gradient vectors for GE scanners
#
# Please use the following configuration on the MR Operator Console:
# - Advanced tab: CV11 - tensor file name: n=4321
# - Diffusion tab: b-value: 3000
# - Diffusion tab: number of T2 : 1 (with 102, 6 x interleaved b0 included in diffusion) or 7 (with 96)
# - Diffusion tab: number of directions (TENSOR): 96 or 102
#
# This file should be placed in /usr/g/bin
#
# Multi-shell details:
# Diffusion 102 (WITH Interleaved b0 values)
# - Interleaved b0: b value = 0 s/mm2; directions = 6
# - Shell 1: b value = 500 s/mm2; directions = 6
# - Shell 2: b value = 1000 s/mm2; directions = 15
# - Shell 3: b value = 2000 s/mm2; directions = 15
# - DTI shell 4: b value = 3000 s/mm2; directions = 60
# Diffusion 96
# - Shell 1: b value = 500 s/mm2; directions = 6
# - Shell 2: b value = 1000 s/mm2; directions = 15
# - Shell 3: b value = 2000 s/mm2; directions = 15
# - Shell 4: b value = 3000 s/mm2; directions = 60
#
# [Reference]
# ABCD Summary of Neuroimaging Parameters
# https://abcdstudy.org/images/Protocol\_Imaging\_Sequences.pdf
```



tensor1127: ADNI-3 Advanced

ADNI (Alzheimer's Disease Neuroimaging Initiative)

tensor1127.dat

- Three-shells DTI
- To match ADNI, use 1 b=0 volume
- # of directions = 126, Max b-value = 2000
 - The “classic” ADNI-3 Advanced
 - 0 x 12
 - 500 x 6
 - 1000 x 48
 - 2000 x 60
 - # of directions = 100, Max b-value = 2000
 - The mayo-modified version of ADNI-3 Advanced for wide bore systems
 - 0 x 9
 - 500 x 6
 - 1000 x 38
 - 2000 x 47



```
# Multi-shell tensor1127.dat file for ADNI-3 Advanced Protocols (Alzheimer's Disease Neuroimaging Initiative)
# Diffusion gradient vectors for GE scanners
#
# Creator: Jaemin Shin Rob Reid, reid.robert@mayo.edu, 2020-04-09
#
# Please use the following configuration on the MR Operator Console:
# - Advanced tab: CV11 - tensor file name: n=1127
# - Diffusion tab: b-value: 2000
# - Diffusion tab: number of T2 : 1
# - Diffusion tab: number of directions (TENSOR): 126 or 100
#
# This file should be placed in /usr/g/bin
#
# Multi-shell details:
# Diffusion 126
# The 126 vector version is the "classic" ADNI3 Advanced 3 shell set.
# It is nominally intended for scanners capable of gradient amplitudes > 70 mT/m.
# - Interleaved b0: b value = 0 s/mm2; directions = 12
# - Shell 1: b value = 500 s/mm2; directions = 6
# - Shell 2: b value = 1000 s/mm2; directions = 48
# - Shell 3: b value = 2000 s/mm2; directions = 60
# Diffusion 100
# The 100 vector version is intended to support multishell scanning in a reasonable time
# on scanners with gradient coils with lower maximum strength.
# - Interleaved b0: b value = 0 s/mm2; directions = 9
# - Shell 1: b value = 500 s/mm2; directions = 6
# - Shell 2: b value = 1000 s/mm2; directions = 38
# - Shell 3: b value = 2000 s/mm2; directions = 47
#
# Both versions have their b values interspersed through the whole scan
# to minimize gradient heating, and the directions are spread over
# entire spheres to aid eddy current correction.
#
# The directions were optimized using an electrostatic repulsion scheme:
# Caruyer, Emmanuel, Christophe Lenglet, Guillermo Sapiro, and Rachid Deriche.
# "Design of multishell sampling schemes with uniform coverage in diffusion MRI."
# Magnetic Resonance in Medicine 69, no. 6 (2013): 1534-1540.
#
# http://adni.loni.usc.edu/adni-3/
#
```



tensor521: UK Biobank

UK Biobank Brain Imaging

tensor521.dat

- Two-shells DTI
- To match UK Biobank, use 1 b=0 volume
- # of directions = 104, Max b-value = 2000
 - Primary directions
 - 0 x 4
 - 1000 x 50
 - 2000 x 50
- # of directions = 6, Max b-value = 2000
 - Secondary directions to get blip-reversed b=0 scans for distortion correction
 - 0 x 3
 - 2000 x 3



```
# Multi-shell tensor521.dat file for UK Biobank Brain Imaging
# Diffusion gradient vectors for GE scanners
#
# Creator: Jaemin Shin (Jaemin.Shin@ge.com), 2020-05-04
#
# Please use the following configuration on the MR Operator Console:
# - Advanced tab: CV11 - tensor file name: n=521
# - Diffusion tab: b-value: 2000
# - Diffusion tab: number of T2 : 1 (4 x interleaved b0 included in diffusion)
# - Diffusion tab: number of directions (TENSOR): 104
#
# This file should be placed in /usr/g/bin
#
# Multi-shell details:
# Diffusion 104
# - Interleaved b0: b value = 0 s/mm2; directions = 4
# - Shell 1: b value = 1000 s/mm2; directions = 50
# - Shell 2: b value = 2000 s/mm2; directions = 50
# Diffusion 6 (optional)
# Short acquisition to get blip-reversed b=0 scans for distortion correction.
# - Interleaved b0: b value = 0 s/mm2; directions = 1
# - Shell 1: b value = 2000 s/mm2; directions = 3
#
# [Reference]
# UK Biobank Brain Imaging - Online Resources
# https://www.fmrib.ox.ac.uk/ukbiobank/protocol/
#
# Diffusion gradient vectors from Siemens DiffusionVectors:
# Primary directions 104 directions
# https://www.fmrib.ox.ac.uk/ukbiobank/protocol/NewDiffusionVectors230914.txt
# Secondary directions 6 directions
# https://www.fmrib.ox.ac.uk/ukbiobank/protocol/DiffusionVectors3b0.txt
#
```



```
# Multi-shell tensor1225.dat file for HCP Lifespan
# Diffusion gradient vectors for GE scanners
# from Lifespan HCP-Development and HCP-Aging protocol package
# (Version 2019.01.14)
#
# 01/25/2022
# Creator: Jaemin Shin(Jaemin.Shin@ge.com), Baolian Yang(Baolian.Yang@ge.com)
#
# Please use the following configuration on the MR Operator Console:
# - Advanced tab: CV11 - tensor file name: n=1225
# - Diffusion tab: b-value: 3000
# - Diffusion tab: number of directions (TENSOR): 198
#
# Multi-shell details:
# Diffusion 198 = Diffusion 98 + 1 x b0 + Diffusion 99
# - Interspersed b0: b value = 0 s/mm2; directions = 13
# - Shell 1: b value = 1500 s/mm2; directions = 92
# - Shell 2: b value = 3000 s/mm2; directions = 93
# Diffusion 98
# - Interspersed b0: b value = 0 s/mm2; directions = 6
# - Shell 1: b value = 1500 s/mm2; directions = 46
# - Shell 2: b value = 3000 s/mm2; directions = 46
# Diffusion 99
# - Interspersed b0: b value = 0 s/mm2; directions = 6
# - Shell 1: b value = 1500 s/mm2; directions = 47
# - Shell 2: b value = 3000 s/mm2; directions = 46
#
```

```
# [Reference]
# Harms MP, Somerville LH, Ances BM, Andersson J, Barch DM, Bastiani M,
# Bookheimer SY, Brown TB, Buckner RL, Burgess GC, Coalson TS.
# Extending the Human Connectome Project across ages: Imaging protocols
# for the Lifespan Development and Aging projects.
# Neuroimage. 2018 Dec 1;183:972-84.
# https://doi.org/10.1016/j.neuroimage.2018.09.060
#
# Emmanuel Caruyer, Christophe Lenglet, Guillermo Sapiro, Rachid Deriche.
# Design of multishell sampling schemes with uniform coverage in diffusion MRI.
# Magnetic Resonance in Medicine, 2013, 69:1534-1540.
# http://dx.doi.org/10.1002/mrm.24736
#
# Tool available at:
# http://www.emmanuelcaruyer.com/q-space-sampling.php
#
# Lifespan HCP-Development and HCP-Aging protocol package at:
# https://www.humanconnectome.org/study/hcp-lifespan-development/project-protocol/imaging-protocols-hcp-development
#
# https://www.humanconnectome.org/storage/app/media/protocols/HCP\_VE11C\_Prisma\_2019.01.14.zip
```

